



To whom it may concern

June 4th , 2024

Temperature during transport of FertiPro products

The following products have been proven suitable for transport or short-term storage at elevated temperatures (up to 5 days at 37°C):

- FertiCult Flushing medium
- Sil-Select Plus Washing/Insemination medium
- FertiCult IVF medium
- FertiCult Aspiration medium
- GAIN medium
- FertiCult Mineral Oil
- FertiCult High Viscosity Oil
- VitriFreeze/VitriThaw
- VitriFreeze ES/VitriThaw ES
- FertiVit Cooling/FertiVit Warming
- SpermFreeze
- SpermFreeze SSP
- Hyaluronidase in FertiCult Flushing medium
- 10% PVP in FertiCult Flushing medium
- Sil-Select Plus
- Sil-Select Stock
- EpiScreen Plus
- SpermMar Test IgG
- SpermMar Test IgA
- LeucoScreen Plus
- VitalScreen
- HOS Test
- Spermac Stain

The following products have been proven suitable for transport or short-term storage at elevated temperatures (up to 5 days at 25°C):

- Fructose Test

All documentation and records of test data have been reviewed by the relevant authorized persons and are retained as specified in ISO 13485.

Signature:

A handwritten signature in blue ink, appearing to read "Ilse Jonckheere".

Ilse Jonckheere
Vice President
FertiPro NV

Для предъявления по месту требования

Июнь 4, 2024

Температурный режим во время транспортировки продукции FertiPro

Следующие продукты сохраняют стабильность при транспортировке и кратковременном хранении при колебаниях температуры до 5 дней до 37 C°:

- FertiCult Flushing medium
- Sil-Select Plus Washing/Insemination medium
- FertiCult IVF medium
- FertiCult Aspiration medium
- GAIN medium
- FertiCult Mineral Oil
- FertiCult High Viscosity Oil
- VitriFreeze/Thaw
- VitriFreeze/Thaw ES
- FertiVit Cooling/Warming
- EmbryoFreeze/Thaw
- SpermFreeze
- SpermFreeze SSP
- Hyaluronidase in FertiCult Flushing medium
- 10% PVP in FertiCult Flushing medium
- Sil-Select Plus
- Sil-Select Stock
- EpiScreen Plus
- SpermMar Test IgG
- SpermMar Test IgA
- LeucoScreen Plus
- VitalScreen
- HOS test
- Spermac Stain

Следующие продукты сохраняют стабильность при транспортировке и кратковременном хранении при колебаниях температуры до 5 дней до 25 C°:

- Fructose Test

Вся документация и отчеты испытаний были рассмотрены соответствующими уполномоченными лицами и соответствуют ISO 13485.

Подпись:

Вице президент FertiPro

Илзе Джонкхир